

Spencer AMBUJET Portable Suction Device

Rechargeable Battery Suction Unit With Wall Mount



Product Overview

Space-saving suction device to be used to remove fluids from the airway or respiratory support system. The device creates a negative pressure that draws fluids through disposable tubing that is connected to the collection container. The base performs the dual function of charging station and fixing in accordance to EN 1789 standard. Securely fitted in a 10G or 20G tested wall mount, this battery-powered suction unit can be removed from the wall with the push of a button.

Key Features

- Low weight and size
- Easy view of the setted vacuum thanks to the led indicators
- Classification ISO High Vacuum/ High Flow
- Suction device protected from bumps and falls thanks to the protective shell
- Mechanism of engagement and release from the base extremely rapid
- Recharging of the battery is done simply by placing the device on the mounting system

Applications

- EMS
- Rescue



www.frsa.com.au | admin@frsa.com.au | + 61 7 3209 7422

Spencer AMBUJET

Portable Suction Device

Rechargeable Battery Suction Unit With Wall Mount



Technical Data	
Jet Dimensions	200x225x285 mm
Jet Dimensions w/fixing	200x250x320 mm
Suction device weight	2.5 ± 0.1 kg
Fixing weight	0.8 ± 0.1kg
Vacuum Range	From 50 to 550 mmHg
Flow (free flow)	32 l/min
Noise level	70 db
Electrical connection requirements	100-240 VAC 47-63 Hz 0,75A max; 12VDC 33W max
IP Grade	IP12
Battery	12V DC Ni-MH
Battery operating time	30 min. ON, 30 min. OFF
Battery recharging time	5 h
Container capacity	1000 ml
Patient tube length	1,8 m
Operating temp.	From 0°C to +45°C
Storage temp.	From – 40°C to +70°C
Humidity	From 0 to 95%
Compliant with	ISO 10079-1 EN 1789
Class IIa Med.Device according to Dir. 93/42/CEE	AMBUJET P 1000 B 20G SC75510D

Standard Equipment
Fixing with integrated charging system
Battery
Serres canister with disposable bag
Disposable patient tube 1,8 m
Shoulder strap for easy transport

Optional Accessories / Consumables	
SC75012E	Lateral bags for Jet
SC73017E	220V Power supply
SC73013E	Antibacterial filter
SC73016E	Sterile suction catheter 1,8m
SC75018E	36 disposable bags for canister



www.frsa.com.au | admin@frsa.com.au | + 61 7 3209 7422